

Product Certificates



Item no.	Description	Cleanroom Class 7	Irradiation	FDA 21CFR §	Food Contact EU 1935/2004 EU 10/2011 raw material/ tested	United States Pharmacopeia USP & European Pharmacopoeia EP	Bisphenol A added/ tested	Phthalates added/ tested	Latex added/ tested	BSE/TSE material added/ EMEA/410/01 fulfilled	Reach (EG) Nr. 1907/2006	RoHS EU 2011/65 conformity of base material	Material
5378-0001	Scoop LaboPlast, single-use, PS, 25 ml, without cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0002	Scoop LaboPlast, single-use, PS, 25 ml, with cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0003	Scoop LaboPlast, single-use, PS, 50 ml, without cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0004	Scoop LaboPlast, single-use, PS, 50 ml, with cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0005	Scoop LaboPlast, single-use, PS, 100 ml, without cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0006	Scoop LaboPlast, single-use, PS, 100 ml, with cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0007	Scoop LaboPlast, single-use, PS, 150 ml, without cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0008	Scoop LaboPlast, single-use, PS, 150 ml, with cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0009	Spatula LaboPlast, single-use, PS, 150 mm, without cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0010	Spatula LaboPlast, single-use, PS, 150 mm, with cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0011	Spoon LaboPlast, single-use, PS, 2,5 ml	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0012	Spoon LaboPlast, single-use, PS, 10 ml	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0013	Scoop LaboPlast, single-use, PS, 250 ml, without cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0014	Scoop LaboPlast, single-use, PS, 250ml, with cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0015	Scoop LaboPlast, single-use, PS, 500 ml, without cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0016	Scoop LaboPlast, single-use, PS, 500ml, with cover	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0019	Scoop LaboPlast, single-use, PS, 1000 ml	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0030	Scraper LaboPlast, single-use, PS, 80 mm	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0031	Spoon spatula, LaboPlast, PS, 0,5 ml/17 mm	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0032	Micro spatula, LaboPlast, PS, 1,0/0,25 ml	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0040	Forceps, disposable, pointed/angled, 130mm, PS	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0041	Forceps, disposable, pointed, 130mm, PS	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-0042	Forceps, disposable, broad tip, 130mm, PS	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1001	Sampling scoop SteriPlast, PS, sterile, 25 ml, without cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1002	Sampling scoop SteriPlast, PS, sterile, 25 ml, with cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1003	Sampling scoop SteriPlast, PS, sterile, 50 ml, without cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1004	Sampling scoop SteriPlast, PS, sterile, 50 ml, with cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1005	Sampling scoop SteriPlast, PS, sterile, 100 ml, without cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1006	Sampling scoop SteriPlast, PS, sterile, 100 ml, with cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1007	Sampling scoop SteriPlast, PS, sterile, 150 ml, without cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1008	Sampling scoop SteriPlast, PS, sterile, 150ml, with cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1009	Sampling spatula SteriPlast, PS, sterile, 150 mm, without cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1010	Sampling spatula SteriPlast, PS, sterile, 150 mm, with cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1011	Sampling spoon SteriPlast, PS, sterile, 2,5 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1012	Sampling spoon SteriPlast, PS, sterile, 10 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1013	Sampling scoop SteriPlast, PS, sterile, 250 ml, without cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1014	Sampling scoop SteriPlast, PS, sterile, 250 ml, with cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1015	Sampling scoop SteriPlast, PS, sterile, 500 ml, without cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1016	Sampling scoop SteriPlast, PS, sterile, 500 ml, with cover	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1019	Sampling scoop SteriPlast, PS, sterile, 1000 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1030	Scraper SteriPlast, PS, sterile, 80 mm	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1031	Spoon spatula, SteriPlast, PS, sterile, 0,5ml/17mm	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1032	Micro-spatula, SteriPlast, PS, sterile, 1,0/0,25ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1040	Forceps, disposable, pointed/angled, PS, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1041	Forceps, disposable, pointed, 130mm, PS, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1042	Forceps, disposable, broad tip, 130mm, PS, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1050	Single-use tweezers, pointed/angled, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1051	Single-use tweezers, pointed, PS, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1052	Single-use tweezers, broad tip, PS, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1060	Detectable forceps, pointed/angled, PS, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1061	Detectable forceps, pointed, PS, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-1062	Detectable forceps, broad tip, PS, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-2001	Detectable scoops, blue, PS, sterile, 25 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-2003	Detectable scoops, blue, PS, sterile, 50 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-2005	Detectable scoops, blue, PS, sterile, 100 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-2007	Detectable scoops, blue, PS, sterile, 150 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-2009	Detectable spatuals, blue, PS, sterile, 150 mm	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-2011	Detectable spoons, blue, PS, sterile, 2,5 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS

Product Certificates



Item no.	Description	Cleanroom Class 7	Irradiation	FDA 21CFR §	Food Contact EU 1935/2004 EU 10/2011 raw material/ tested	United States Pharmacopeia USP & European Pharmacopoeia EP	Bisphenol A added/ tested	Phthalates added/ tested	Latex added/ tested	BSE/TSE material added/ EMEA/410/01 fulfilled	Reach (EG) Nr. 1907/2006	RoHS EU 2011/65 conformity of base material	Material
5378-2012	Detectable spoons, blue, PS, sterile, 10 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-2013	Detectable scoops, blue, PS, sterile, 250 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-2015	Detectable scoops, blue, PS, sterile, 500 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-2019	Detectable scoops, blue, PS, sterile, 1000 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-2030	Detectable scrapers, blue, PS, sterile, 80 mm	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3001	Food scoops, blue, PS, sterile, 25 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3003	Food scoops, blue, PS, sterile, 50 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3005	Food scoops, blue, PS, sterile, 100 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3007	Food scoops, blue, PS, sterile, 150 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3009	Food spatulas, blue, PS, sterile, 150 mm	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3011	Food spoons, blue, PS, sterile, 2,5 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3012	Food spoons, blue, PS, sterile 10 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3013	Food scoops, blue, PS, sterile, 250 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3015	Food scoops, blue, PS, sterile, 500 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3019	Food scoops, blue, PS, sterile, 1000 ml	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-3030	Food scrapers, blue, PS, sterile, 80 mm	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4004	Dosing spoons 0.5 ml, PS, white	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4005	Dosing spoons 0.5 ml, PS, white, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4006	Dosing spoons 1.0 ml, PS, white	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4007	Dosing spoons 1.0 ml, PS, white, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4008	Dosing spoons 2.5 ml, PS, white	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4009	Dosing spoons 2.5 ml, PS, white, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4010	Dosing spoons 5.0 ml, PS, white	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4011	Dosing spoons 5.0 ml, PS, white, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4012	Dosing spoons 10 ml, PS, white	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4013	Dosing spoons 10 ml, PS, white, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4014	Dosing spoons 15 ml, PS, white	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4015	Dosing spoons 15 ml, PS, white, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4016	Dosing spoons 25 ml, PS, white	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4017	Dosing spoons 25 ml, PS, white, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4018	Dosing spoons 50 ml, PS, white	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4019	Dosing spoons 50 ml, PS, white, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4044	Measuring spoons 0.5 ml, PS, blue	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4045	Measuring spoons 0.5 ml, PS, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4046	Measuring spoons 1.0 ml, PS, blue	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4047	Measuring spoons 1.0 ml, PS, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4048	Measuring spoons 2.5 ml, PS, blue	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4049	Measuring spoons 2.5 ml, PS, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4050	Measuring spoons 5.0 ml, PS, blue	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4051	Measuring spoons 5.0 ml, PS, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4052	Measuring spoons 10 ml, PS, blue	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4053	Measuring spoons 10 ml, PS, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4054	Measuring spoons 15 ml, PS, blue	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4055	Measuring spoons 15 ml, PS, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4056	Measuring spoons 25 ml, PS, blue	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4057	Measuring spoons 25 ml, PS, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4058	Measuring spoons 50 ml, PS, blue	yes	no	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4059	Measuring spoons 50 ml, PS, blue, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4085	Volumetric spoon 0,5ml, PS blue detectable,sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4087	Volumetric spoon 1,0ml, PS blue detectable,sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4089	Volumetric spoon 2,5ml, PS blue detectable,sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4091	Volumetric spoon 5,0ml, PS blue detectable,sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4093	Volumetric spoon 10ml, PS blue detectable, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4095	Volumetric spoon 15ml, PS blue detectable, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4097	Volumetric spoon 25ml, PS blue detectable, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5378-4099	Volumetric spoon 50ml, PS blue detectable, sterile	yes	gamma	§177.1640	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Polystyrene PS
5379-0001	Scoop LaboPlast Bio, single-use, Bio-PE, 25 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-0003	Scoop LaboPlast Bio, single-use, Bio-PE, 50 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE

Product Certificates



Item no.	Description	Cleanroom Class 7	Irradiation	FDA 21CFR §	Food Contact EU 1935/2004 EU 10/2011 raw material/ tested	United States Pharmacopeia USP & European Pharmacopoeia EP	Bisphenol A added/ tested	Phthalates added/ tested	Latex added/ tested	BSE/TSE material added/ EMEA/410/01 fulfilled	Reach (EG) Nr. 1907/2006	RoHS EU 2011/65 conformity of base material	Material
5379-0005	Scoop LaboPlast Bio, single-use, Bio-PE, 100 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-0007	Scoop LaboPlast Bio, single-use, Bio-PE, 150 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-0009	Spatula LaboPlast Bio, single-use, Bio-PE, 150 mm	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-0011	Spoon LaboPlast Bio, single-use, Bio-PE, 2,5 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-0012	Spoon LaboPlast Bio, single-use, Bio-PE, 10 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-1001	Sampling scoop SteriPlast Bio, Bio-PE, sterile, 25 ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-1003	Sampling scoop SteriPlast Bio, Bio-PE, sterile, 50 ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-1005	Sampling scoop SteriPlast Bio, Bio-PE, sterile, 100 ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-1007	Sampling scoop SteriPlast Bio, Bio-PE, sterile, 150 ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-1009	Sampling spatula SteriPlast Bio, Bio-PE, sterile, 150 mm	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-1011	Sampling spoon SteriPlast Bio, Bio-PE, sterile, 2,5 ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5379-1012	Sampling spoon SteriPlast Bio, Bio-PE, sterile, 10 ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	Bio-HDPE
5393-1120	LiquiDispo LaboPlast, HDPE, 500 mm, 100 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-1121	LiquiDispo SteriPlast, HDPE, sterile, 500mm, 100ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-1130	LiquiDispo LaboPlast, HDPE, 1000 mm, 200 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-1131	LiquiDispo SteriPlast, HDPE sterile, 1000mm, 200ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-2240	MicroDispo LaboPlast, HDPE, 500 mm, 10 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-2241	MicroDispo SteriPlast, HDPE, sterile, 500 mm, 10ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-2250	MicroDispo LaboPlast, HDPE, 1000 mm, 10 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-2251	MicroDispo SteriPlast, HDPE, sterile, 1000mm, 10ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-3320	ViscoDispo LaboPlast, HDPE, 500 mm, 100 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-3321	ViscoDispo SteriPlast, HDPE, sterile, 500mm, 100ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-3330	ViscoDispo LaboPlast, HDPE, 1000 mm, 200 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-3331	ViscoDispo SteriPlast, HDPE sterile, 1000mm, 200ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4420	ZoneDispo LaboPlast, HDPE, 500 mm, 100 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4421	ZoneDispo SteriPlast, HDPE, sterile, 500 mm, 100ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4430	ZoneDispo LaboPlast, HDPE, 1000 mm, 200 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4431	ZoneDispo SteriPlast, HDPE, sterile, 1000mm, 200ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4440	MultiDispo LaboPlast, HDPE, 500 mm, 100 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4441	MultiDispo SteriPlast, HDPE, sterile, 500mm, 100ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4450	MultiDispo LaboPlast, HDPE, 1000 mm, 200 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4451	MultiDispo SteriPlast, HDPE sterile, 1000mm, 200ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4460	TargetDispo LaboPlast, HDPE, 500 mm, 100 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4461	TargetDispo SteriPlast, HDPE sterile, 500mm, 100ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4470	TargetDispo LaboPlast, HDPE, 1000 mm, 200 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-4471	TargetDispo SteriPlast, HDPE sterile, 1000mm,200ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5500	DispoPipette LaboPlast, HDPE, 500 mm, 100 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5501	DispoPipette SteriPlast, HDPE sterile, 500mm,100ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5510	PowderDispo LaboPlast, HDPE, 300 mm, 60 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5511	PowderDispo SteriPlast, HDPE, sterile, 300mm, 60ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5521	PowderDispo SteriPlast, HDPE sterile, 500mm, 100ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5522	PowderDispo LaboPlast, HDPE, 500 mm, 100 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5531	DispoPipette LaboPlast, HDPE, 1000 mm, 200 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5532	DispoPipette SteriPlast, HDPE sterile,1000mm,200ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5540	DispoPicker LaboPlast, HDPE, 300 mm, 60 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5541	DispoPicker SteriPlast, HDPE, sterile, 300mm, 60ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5550	DispoPicker LaboPlast, HDPE, 500 mm, 100 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5551	DispoPicker SteriPlast, HDPE sterile, 500mm, 100ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5560	DispoLance LaboPlast, HDPE, 300 mm	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5561	DispoLance SteriPlast, HDPE, sterile, 300 mm	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5570	DispoLance LaboPlast, HDPE, 500 mm	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-5571	DispoLance SteriPlast, HDPE, sterile, 500 mm	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-6620	DispoTube LaboPlast, HDPE, 500 mm, 100 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-6621	DispoTube SteriPlast, HDPE, sterile, 500 mm, 100ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-6630	DispoTube LaboPlast, HDPE, 1000 mm, 200 ml	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-6631	DispoTube SteriPlast, HDPE, sterile, 1000mm, 200ml	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-6640	DispoLadle LaboPlast, HDPE, 500 mm	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-6641	DispoLadle SteriPlast, HDPE, sterile, 500 mm	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE

Product Certificates



Item no.	Description	Cleanroom Class 7	Irradiation	FDA 21CFR §	Food Contact EU 1935/2004 EU 10/2011 raw material/ tested	United States Pharmacopeia USP & European Pharmacopoeia EP	Bisphenol A added/ tested	Phthalates added/ tested	Latex added/ tested	BSE/TSE material added/ EMEA/410/01 fulfilled	Reach (EG) Nr. 1907/2006	RoHS EU 2011/65 conformity of base material	Material
5393-6650	DispoLadle LaboPlast, HDPE, 1000 mm	yes	no	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE
5393-6651	DispoLadle SteriPlast, HDPE, sterile, 1000 mm	yes	gamma	§177.1520	yes / yes	not tested	no/no	no/no	no/no	no/yes	yes	yes	HDPE

The technical information provided herein is presented in good faith, to the best of our knowledge and based on information we believe it is reliable.

The customer is responsible for the appropriate, safe and legal use, processing and handling of our products.

Haltbarkeitsdatum/Shelf life/Péremption: Production date (Lot No.: YYMMDD____xxxxxxx) + 2 years - (Please note: This only refers to the sterility of the SteriPlast products)

This document was issued electronically and is therefore valid without signature.